



COVERS WHEELCHAIR RELATED



POLYMAILLE®

- Coating of polyurethane on Jersey polyester
- Integral cover over 6 sides
- Non-slip base
- Adapted for home care and medical institutions use
- Hypoallergenic
- Anti-bacterial, anti-germ, anti-fungal
- bi-extensible
- Hyper-breathable (liquid repellent material and permeable to steam), adapted for incontinence. Reduce maceration effects through easing the skin respiration.
- Zip on 3 sides (except P311C)
- Integrated handle




POLYMAILLE® HD

- Coating of polyurethane on Jersey polyamide
- Integral cover over 6 sides
- Adapted for heavy care and medical institutions use
- Hypoallergenic
- Anti-bacterial, anti-germ, anti-fungal
- bi-extensible
- Hyper-breathable (liquid repellent material and permeable to steam), adapted for incontinence. Reduce maceration effects through easing the skin respiration.
- Welded and sewed
- Zip on 2 sides (except P311C)
- With screen printing

	TYPE DE PROTECTION	
	POLYMAILLE®	POLYMAILLE® HD
Composition	Coating of Polyurethane on Jersey polyester	Coating of Polyurethane on Jersey polyamide with printed instructions
Version	Integral cover	Integral cover
MATERIAL CHARACTERISTICS		
Anti-allergenic	✓	✓
Antibacterial	✓	✓
Antimicrobial	✓	✓
Antifungal	✓	✓
Bi-extensibility	+	++
Respirability	++	++
TAILORING CHARACTERISTICS (Cushions only)		
Kind of assembly	Sewing	Welded seams
Anti-slip basis	✓	
Opening	3 sides (Except P311C: lower face)	2 sides
Integrated handle	✓	
Serigraphy		Printed washing and disinfection instructions
MAINTENANCE & DISINFECTION		
Wash	95° 1000 RPM	95° 1000 RPM
Disinfection	Disinfect using cold-application spray products (Quaternary ammonium - glutaraldéhydes...)	

GOOD TO KNOW... ON THE COVER'S LABEL



REPLACEMENT COVERS WHEELCHAIR RELATED

DESIGNATION	ITEM CODE POLYMAILLE®	ITEM CODE POLYMAILLE®HD	DIMENSIONS
Cover for SYST'AM® P361CLX Visco-Foam Cushion	P514PA4040W		40 x 40 x 6 cm
	P514PA4242W		42 x 42 x 6 cm
	P514PA4345W		43 x 45 x 6 cm
	P514PA4348W		43 x 48 x 6 cm
	P514PA4543W		45 x 43 x 6 cm
	P514PA4843W		48 x 43 x 6 cm
Cover for SYST'AM® P311C / Gel or P311C-D / SYST'AM® DUOGEL®	P511PA4040W	P511PA4040F	40 x 40 x 2,5 cm
	P511PA4543W	P511PA4543F	45 x 43 x 2,5 cm
Cover for SYST'AM® P313C / Gel with removable central pad	P513PA4040W		40 x 40 x 2,5 cm
Cover for SYST'AM® P314C / MIXTE Cushion	P514PA4040W	P514PA4040F	40 x 40 x 6 cm
	P514PA4543W	P514PA4543F	45 x 43 x 6 cm
Cover for SYST'AM® P341C / DUOFORM® or SYST'AM® P371C / HR-Foam Cushion	P541PA3636W		36 x 36 x 8 cm
	P541PA4242W	P541PA4242F	42 x 42 x 8 cm
	P541PA4246W		42 x 46 x 8 cm
	P541PA4542W	P541PA4542F	45 x 42 x 8 cm
	P541PA4546W	P541PA4546F	45 x 46 x 8 cm
	P541PA5146W	P541PA5146F	51 x 46 x 8 cm
Cover for SYST'AM® P361C / VISCOFLEX®	P541PA5150W		51 x 50 x 8 cm
	P541PA3636W		36 x 36 x 8 cm
	P541PA3838W		38 x 38 x 8 cm
	P541PA4040W		40 x 40 x 8 cm
	P541PA4242W		42 x 42 x 8 cm
	P541PA4246W		42 x 46 x 8 cm
	P541PA4542W		45 x 42 x 8 cm
	P541PA4546W		45 x 46 x 8 cm
	P541PA4843W		48 x 43 x 8 cm
Cover for SYST'AM® P361CA / VISCOFLEX®+	P541PA5146W		51 x 46 x 8 cm
	P541PA5150W		51 x 50 x 8 cm
	P541PA3232W		32 x 32 x 8 cm
	P541PA3636W		36 x 36 x 8 cm
	P541PA3838W		38 x 38 x 8 cm
	P541PA3843W		38 x 43 x 8 cm
	P541PA4040W	P541PA4040F	40 x 40 x 8 cm
	P541PA4242W	P541PA4242F	42 x 42 x 8 cm
	P541PA4246W		42 x 46 x 8 cm
	P541PA4250W		42 x 50 x 8 cm
	P541PA4542W	P541PA4542F	45 x 42 x 8 cm
	P541PA4546W	P541PA4546F	45 x 46 x 8 cm
	P541PA4550W		45 x 50 x 8 cm
	P541PA4843W	P541PA4843F	48 x 43 x 8 cm
Cover for SYST'AM® P361CA / VISCOFLEX®+ Bariatric	P541PA5146W	P541PA5146F	51 x 46 x 8 cm
	P541PA5150W	P541PA5150F	51 x 50 x 8 cm
Cover for SYST'AM® P361CP / VISCOFLEX®+ with pommel	P541PA5646W	P541PA5646F	56 x 46 x 10 cm
	P541PA6150W		61 x 50 x 10 cm
	P541PAP4242W		42 x 42 x 8 cm
	P541PAP4542W		45 x 42 x 8 cm
	P541PAP4546W		45 x 46 x 8 cm

DESIGNATION	ITEM CODE POLYMAILLE®	ITEM CODE POLYMAILLE®HD	DIMENSIONS
Cover for SYST'AM® P361CR / VISCOFLEX® Air	P502PA4343W	P502PA4343F	43 x 43 x 9,5 cm
	P502PA4643W	P502PA4643F	46 x 4 x 9,5 cm
	P502PA4646W	P502PA4646F	46 x 46 x 9,5 cm
Cover for SYST'AM® P301C / POLYAIR® 60	P501PA3636W		36 x 36 x 6 cm
	P501PA4040W		40 x 40 x 6 cm
	P501PA4043W		40 x 43 x 6 cm
	P501PA4046W		40 x 46 x 6 cm
	P501PA4343W		43 x 43 x 6 cm
	P501PA4346W		43 x 46 x 6 cm
	P501PA4640W		46 x 40 x 6 cm
	P501PA4646W		46 x 46 x 6 cm
	P501PA5146W		51 x 46 x 6 cm
	P501PA5150W		51 x 50 x 6 cm
	P501PA5646W		56 x 46 x 6 cm
	P501PA6150W		61 x 50 x 6 cm
Cover for SYST'AM® P302C / POLYAIR® 100	P502PA4040W		40 x 40 x 10 cm
	P502PA4043W		40 x 43 x 10 cm
	P502PA4046W		40 x 46 x 10 cm
	P502PA4343W		43 x 43 x 10 cm
	P502PA4346W		43 x 46 x 10 cm
	P502PA4640W		46 x 40 x 10 cm
	P502PA4646W		46 x 46 x 10 cm
	P502PA5146W		51 x 46 x 10 cm
	P502PA5150W		51 x 50 x 10 cm
	P502PA5646W		56 x 46 x 10 cm
Cover for SYST'AM® P951D / BACK SUPPORT	Small	P951SHW	44 x 39 x 9 cm
	Standard	P951HW	49 x 44 x 9 cm
	Large	P951LHW	49 x 50 x 10 cm
Cover for SYST'AM® P961F / UPPER LIMBS POSITIONING WEDGE IN CHAIR		P961HW	

PACKAGING



		DIMENSIONS
REGULAR CUSHIONS SIZES		
32 x 32 cm, 36 x 36 cm, 38 x 38 cm, 38 x 43 cm, 40 x 40 cm, 40 x 43 cm, 40 x 46 cm, 42 x 42 cm, 42 x 46 cm, 43 x 43 cm, 43 x 46 cm, 45 x 42 cm, 45 x 43 cm, 45 x 46 cm, 46 x 40 cm, 46 x 46 cm, 48 x 43 cm	/1 Standard box	48 x 43,5 x 10,5 cm
	/4	53 x 47 x 42,5 cm
LARGE CUSHIONS SIZES AND OTHER ITEMS		
42 x 50 cm, 45 x 50 cm, 51 x 46 cm, 51 x 50 cm, 56 x 46 cm, 61 x 50 cm	/1 XL Box	63 x 54,5 x 10,5 cm
	/2	67,5 x 56 x 21 cm
P951D Back support	/1	Small : 48 x 43,5 x 10,5 cm Standard : 48 x 43,5 x 10,5 cm Large : 63 x 54,5 x 10,5 cm
P952D Back support	/1	53 x 47 x 11,5 cm
P961F Upper limbs positioning wedge in chair	/1	23 x 19 x 52 cm